

P91350Hu02
Adipose Differentiation Related Protein (ADRP)
Organism: Homo sapiens (Human)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[DESCRIPTION]

Protein Names: Adipose Differentiation Related Protein

Synonyms: ADRP, PLIN2,ADFP

Species: Human

Size: 100µg

Source: *Escherichia coli*-derived

Subcellular Location: Membrane; Peripheral membrane protein.

[PROPERTIES]

Residues: Lys245~Glu433 (Accession # Q99541), with N-terminal His-Tag.

Grade & Purity: >95%, 24kDa as determined by SDS-PAGE reducing conditions.

Formulation: Supplied as lyophilized form in PBS, pH 7.4, containing 5% sucrose, 0.01% sarcosyl.


Endotoxin Level: <1.0 EU per 1µg (determined by the LAL method).

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 23.0kDa

Predicted isoelectric point: 5.8





[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MGHHHHHSGSEF- KSQQTI SQLHSTVHLI EFARKNVYSA NQKIQDAQDK LYLSWVEWKR
SIGYDDTDES HCAEHIESRT LAIARNLTQQ LQTTCHTLLS NIQGVQNIQ DQAKHMGVMA
GDIYSVFRNA ASFKEVSDSL LTSSKGQLQK MKESLDDVMD YLVNNTPLNW LVGPFYPLT
ESQNAQDQGA EMDKSSQETQ RSE

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, BELGIUM
Tel 0032 16 58 90 45 | Fax 0032 16 50 90 45
www.gentaur-worldwide.com
info@gentaur.com