



Affinity Purified Goat anti-Guinea Pig IgG2

Host Goat

Catalog Number GG2-60A

Quantity 1.0 mg

Physical State Liquid

Buffer PBS, pH 7.2

Preservative 0.1% Sodium Azide

Storage 4 C

Procedure

Goats were immunized with highly purified guinea pig IgG2 in Freunds adjuvant. The resulting antiserum was sequentially absorbed with an IgG1 immunoabsorbants until it was IgG2 specific. The anti-Guinea Pig IgG2 was then immunoaffinity purified and tested for monospecificity.

Specification

The anti IgG1 only reacts with IgG1 by IEP and gel diffusion analysis when tested against each of the guinea pig IgG subclasses at 1.0 mg/ml.

Western Blots 1 ug/ml - 0.1 ug/ml

ELISA coating 5 ug/ml

Application For Research and Manufacturing Only; Not for Use in Diagnostic Applications.