



## DATA SHEET

### Product Description

<b>Product:</b>	Affinity Purified Rabbit Anti Rat uPA
<b>Catalog Number:</b>	ASRUPA-GF-HT
<b>Description:</b>	Polyclonal antibody (host rabbit). Affinity purified by immobilized rat uPA.
<b>Recommended Storage:</b>	-70 C
<b>Shipping Conditions:</b>	Dry ice
<b>Minimum Shelf Life:</b>	5 years from delivery

### Technical Information

<b>Purity:</b>	Affinity purification
<b>Source:</b>	Rabbit serum
<b>Molecular Weight:</b>	160,000
<b>Extinction Coefficient:</b>	1.36
<b>Antigen:</b>	Recombinant full length wild-type mouse protein (Glycosylated).
<b>Target:</b>	Recognises Urokinase (uPA)
<b>Target Molecular Weight:</b>	44,000
<b>Host:</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Applications:</b>	ELISA, WB, IHC-P, IHC-Fr
<b>Cross-reactivity:</b>	Rat

### Sample Characteristics

<b>Lot Number:</b>	ASRUPA-GF-HT-SAMPLE
<b>Aliquot:</b>	1 x 0.1 mg
<b>Concentration:</b>	10.0 mg/ml
<b>Volume:</b>	0.01 ml
<b>Buffer:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>Form:</b>	Frozen liquid

Optimal experimental conditions and working dilutions must be determined by the end user. Sample concentration and volume are for example purposes only and do not necessarily represent the product that will be received. Please contact us for a current lot-specific product datasheet.

For *in vitro* laboratory use only

Gentaur Molecular Products  
Voortstraat 49  
1910 Kampenhout, Belgium