Product #	Description
AMRE776	Agarose 3:1 High Resolution Blend, gelling at 34-38°C, EEO ~0.12

· · modeli						
Product:	CERTIFICATE OF QUALITY AGAROSE 3:1					
Grade:	BIOTECHNOLOGY G	GRADE	CAS Number: 9012-36-6 Synonym(s): High Resolution Agarose Blend Molecular Weight: N/A Molecular Formula: N/A			
Appearance:	Fine white powder		Storage Condition: ROOM TEMPERATURE Risk Statements: N/A			
Mol. Formula:			Safety Statements: N/A			
Mol. Weight:	Not Applicable		Hazard Code: None	UN Number: NONE		
Description:						
Assay:		<u>Min</u>	<u>Max</u>	<u>Units</u>	Disposition	
Appearance		PASS		P/F	PASS	
DNase			NONE	P/F	PASS	
EEO (-mr)			0.12		PASS	
Endonuclease/ligase inhibitory factor			NONE	P/F	PASS	
Function Test		PASS		P/F	PASS	
Gel Strength (1	.5%)	2000		g/cm2	PASS	
Gelling Range ((1.5%)	34	38	degrees C	PASS	
Melting Point (1	.5%)	85	89	degrees C	PASS	
Protease			NONE	P/F	PASS	
Reassay Date		REPOR	Т		PASS	
RNase			NONE	P/F	PASS	
Sulfate			0.10	%	PASS	

Storage: Store at room temperature. Protect from moisture.

For research only. Not for human or diagnostic use