

DATA SHEET

Product	Ficoll 400
Grade	ACS
Product Code	F760014
Formula	N/A

MW 350,000-500,000 CAS# 26873-85-8

Lot No

Test Items	Specifications	Results
Appearance	White to off-white crystalling	ie
Solubility	Water-soluble	
Hydrolysis	pH<3	
Stable	pH>=7.0	

Description: A nonionic synthetic polymer of sucrose.

Storage: R.T.

Caution:

Avoid contact with eyes, skin and clothing.

Wash all areas of contact with copious amounts of water. See Material Safety Data Sheet for additional information.

For laboratory or further manufacturing use only.

Not intended for household use.

Laboratory Control