

Herpesvirus Type 1 (BHV-1_IBR), Goat anti-Bovine Cat. # YV0082

**INFECTIOUS BOVINE
RHINOTRACHEITIS
VIRUS/BOVINE HERPESVIRUS
TYPE 1 (IBR/BHV-1)**
Anti-Viral Antiserum

CATALOG NO.: YV0082

VOLUME: 2 ml

LOT:

EXPIRATION:

VIRUS: Infectious Bovine Rhinotracheitis
Virus/ Bovine Herpesvirus Type 1
(IBR/BHV-1)



DESCRIPTION: Anti-IBR polyclonal antiserum. Goat origin.

QUALITY CONTROL METHOD: Indirect FA using IBR FA 12-well Slides;
Rabbit Anti-Goat IgG FITC Conjugate.

Specific Reaction: 2-4+ positive fluorescence at 1:800.

Other Comments: Endpoint titer at 1:3200.

INTENDED USE: Useful for IFA. Not suitable for cell culture serum neutralization because it contains sodium azide. Does not react with other common bovine viruses.

STORAGE: This antiserum is provided in liquid form and should be stored at 4-8°C. DO NOT FREEZE! If antiserum becomes cloudy, it should be discarded. This antiserum contains 0.09% sodium azide as a preservative.