



(CPT1; CPT1-L; L-CPT1; CPT1A)

Product Name	carnitine palmitoyltransferase 1A (liver)
Price	150 € / 500ul
Primary source	HGNC:2328
Gène type	Protéin Coding
See related	Ensembl:ENSG00000110090; HPRD:02755; MIM:600528
Usage	Western Blot (1/500) , ELISA, immunoprecipitation.
Raised in	Rabbit
Clonality	Polyclonal
Form	Crude serum
Storage instructions	Aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Expiration	20/04/2016
Reference	Guido Bommer- Institut de Duve

Official Symbol CPT IA

Official Full Name carnitine palmitoyltransferase 1A (liver)

Organism Homo Sapiens

Lineage : *Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Primates; Haplorrhini; Catarrhini; Hominidae; Homo*

Also known as CPT1; CPT1-L; L-CPT1; CPT1A

