



RNA Templates							DNA Templates							
Probe Assays		SYBR Green Assays					HRM	Probe Assays		SYBR Green Assays				
SensiMix™ Capillary One-Step Kit	SensiMix™ Probe One-Step Kit	SensiMix™ SYBR & Fluorescein One-Step Kit	SensiMix™ Capillary One-Step Kit	SensiMix™ SYBR No-ROX One-Step Kit	SensiMix™ SYBR Low-ROX One-Step Kit	SensiMix™ SYBR One-Step Kit	SensiMix™ HRM Kit	SensiMix™ Capillary Kit	SensiMix™ II Probe Kit	SensiMix™ SYBR & Fluorescein Kit	SensiMix™ Capillary Kit	SensiMix™ SYBR No-ROX Kit	SensiMix™ SYBR Low-ROX Kit	SensiMix™ SYBR Kit

Applied Biosystems™														
7000, 7300, 7700, 7900, 7900HT, Stepone™, Stepone™ plus		●				●			●					●
7500		●			●				●				●	
Roche														
Lightcycler® 1.0, Lightcycler® 2.0	●			●				●			●			
Lightcycler® 480		●		●			●		●			●		
BIO-RAD														
iCycler® MyiQ™, iQ™5		●	●						●	●				
Opticon™, Opticon™2, Chromo4™, MiniOpticon™, CFX96™, CFX384™		●		●					●			●		
Cepheid®														
SmartCycler™		●		●					●			●		
Corbett														
Rotor-Gene™ 3000		●		●					●			●		
Rotor-Gene™ 6000		●		●			●		●			●		
Eppendorf														
Mastercycler® ep realplex		●		●					●			●		
Techne														
Quanta®		●		●					●			●		
Stratagene														
MX4000P™, MX3000P™, MX3005P™		●			●				●			●		