



DATA SHEET

Product Description

Product:	Active Dog Urokinase
Catalog Number:	DUPA
Description:	Active two-chain canine urokinase (uPA), recombinantly produced in insect cells.
Recommended Storage:	-70 C
Shipping Conditions:	Dry ice
Minimum Shelf Life:	3 years from delivery

Technical Information

Purity:	>95% by SDS-PAGE analysis
Source:	Insect cell culture
Molecular Weight:	50,000
Extinction Coefficient:	1.36

Sample Characteristics

Lot Number:	DUPA-SAMPLE
Aliquot:	1 x 1.0 mg
Concentration:	4.0 mg/ml
Volume:	0.25 ml
Buffer:	0.05M Sodium Acetate; 0.1M NaCl; 1mM EDTA; pH 5.0
Form:	Frozen liquid

Optimal experimental conditions and working dilutions must be determined by the end user. Sample concentration and volume are for example purposes only and do not necessarily represent the product that will be received. Please contact us for a current lot-specific product datasheet.

For *in vitro* laboratory use only

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium