



PRODUCT PROFILE

PRODUCT: Affinity Purified Horse anti-Hepatitis B

CATALOG NUMBER: HHPB-55A

LOT NUMBER: 6

QUANTITY: 1.0 mL @ 1.0 mg/mL

PHYSICAL STATE: Liquid

BUFFER: Phosphate buffered saline, pH 7.2

PRESERVATIVE: 0.1% sodium azide

STORAGE: 4° C

EXPIRATION: 24 months from the date the reagent is shipped.

PROCEDURE: A horse was immunized with highly purified native Ad/Ay antigen. Serum was collected and the anti-Ad/Ay was immunoaffinity purified.

SPECIFICATIONS: The anti-hepatitis b is suitable for ELISA, precipitation, fluorescence and immunohistochemistry techniques.

FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION

CERTIFICATE OF ANALYSIS VERIFICATION:
PRODUCT OF THE UNITED STATES OF AMERICA

DATE: _____

Gentaaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium

